


Razer Tartarus Diagram for Synthesis PACS



- 01 → Refresh
- 02 → Add Montage
- 03 → Montage/All Exam
- 04 → Cine
- 05 → Toolbox
- 06 → Edge Enhance
- 07 → Arrow
- 08 → Annotate
- 09 → Spine Labeling
- 10 → Power Scroll
- 11 → ROI
- 12 → Measure
- 13 → Magnify
- 14 → Pan
- 15 → Page/Page All
- 16 → Open Speech Org
- 17 → Submit Speech Org
- 18 → AI (auto impress)
- 19 → Undo (report only)

Sign & Next



Up in Report

First Line of Impression (middle click)

Down in Report

Change Viewport

Scroll Up Images

W/L Favorites

Scroll Down Images

Mic On/Off